


Hydrash



Hydrodermabrasion Technology

Hydrash _____ Hydrodermabrasion Technology

**The newest way to
way to hydrate and
detoxify and detoxify
skin.**



servicioalcliente@luxomedica.com
www.tentrekus.com

Copyright © 2022 Tentrek Lasers. All rights reserved



Hydrash

TIME
15

MIN MAX

UP
DOWN

TIME LEVEL

PULSE

START/STOP POWER

HYDRASH 5.000 cm

UP
DOWN

MODE SELECT

READY/CLEANING

START/STOP POWER

ULTRASOUND

HYDRASH

TECHNICAL

Hydrash, a premium experience


Hydrash is designed to promote exfoliation and removal of aged dead skin cells through hydrodermabrasion.

Using a combination of oxygen and aqueous solutions at supersonic speeds, accumulated debris is removed from the skin, achieving improved luminosity and texture.

Hydrash represents the most advanced technology to stimulate natural skin restoration, providing instant results, without discomfort or downtime.

Hydrash renews the concept of premium facial cleansing, taking it to the next level.

Hydrash





Great results

Reached with their 3 hand pieces

Microdermabrasion


It has two interchangeable tips with diamond crystals. Its function is to remove dead cells, improving circulation and overall skin tone.



Tip with diamond crystals.


Hydrodermabrasion

This handpiece is responsible for supplying water and at the same time sucking impurities from the skin, hydrating it as it goes, giving it a radiant and toxin-free appearance.



Ultrasound

This piece is used to perform sonophoresis, dilating the pores through a cellular massage and then applying active products such as anti-acne peeling or moisturizers.



After the entire treatment, you will experience a feeling of relaxation. Ultrasound can release pain and unclog pores to treat fat.

Advanced Cleaning

The newest way to hydrate your skin

Hydrash technology offers a unique and innovative skin rejuvenation treatment:

1. Skin cleansing, in which it is prepared for treatment.
2. Peeling, in which the active ingredient is placed on the skin and the pores are opened.
3. Suction pump or hydrodermabrasion, in which the superficial layer of the skin is removed (dead cells), impurities and blackheads are eliminated.
4. Sonophoresis, using ultrasound at 1.2 MHz, the anti-acne, anti-grease or anti-aging peel is introduced as required.

Hydrash accelerates the cell renewal process, resulting in smoother skin.

Hydrash




Clinical Cases

Learn about some of the results obtained with Tensamax:



BLACKHEADS

Courtesy: Tentrek Academy



COMEDONES

Courtesy: Dr. Julieth Martinez



BLACKHEADS

Courtesy: Dr. Johvanna Tamayo



ADVANCED CLEANING

Courtesy: Dr. Maria Paula Gaitan



BLACKHEADS

Courtesy: Dr. Freddy Hernandez



Hydrash

— Technical data sheet

| | |
|--------------------------------|-------------------------------|
| Ultrasound piece | 1.2 MHz |
| Hydrodermabrasion piece | Face (regular), nose (narrow) |
| Microdermabrasion piece | #140 / #170 |
| Suction | Until 600mmHg ($\pm 10\%$) |
| Consumption | 80Watts |
| Reservoir capacity | 1200cc |
| Feeding | 110VAC / 50 / 60Hz |
| Dimensions (al*an*prof) | 1000x350x350mm |
| Weight | 30kg |

Hydrash



Hydrash

Registro sanitario:
2018DM-0018354
Resolución:
No. 2018033043

www.tentrekus.com



@tentreklasers