



**Hydra-SH** \_\_\_\_\_  
Hydrodermabrasion Technology

**The newest way  
to hydrate and  
detoxify  
skin.**

[www.tentrekusa.com](http://www.tentrekusa.com)

[www.reviveregenerative.com](http://www.reviveregenerative.com)

**REVIVE**

POWERED BY  **TENTREK**



Hydrazon

15

ON/OFF

100% 11.00 5.00

ULTRASOUND

HYDRASH

UP  
DOWN

UP  
DOWN

TIME LEVEL

MODE SELECT

PULSE

READY/CLEANING

START/STOP POWER

START/STOP POWER

Hydrazon

# Hydra-sh, a premium experience

Hydra-sh is designed to promote exfoliation and removal of aged dead skin cells through hydrodermabrasion.

Using a combination of oxygen and aqueous solutions at supersonic speeds, accumulated debris is removed from the skin, achieving **improved luminosity** and **texture**.

Hydra-sh represents the most advanced technology to stimulate **natural skin restoration**, providing instant results, without discomfort or downtime.

Hydra-sh renews the concept of premium facial cleansing, taking it to the next level.

**Hydra-SH**





# Great results

Reached with their 3 hand pieces

## Microdermabrasion


It has two interchangeable tips with diamond crystals. Its function is to remove dead cells, improving circulation and overall skin tone.



Tip with diamond crystals.


## Hydrodermabrasion

This handpiece is responsible for supplying water and at the same time sucking impurities from the skin, hydrating it as it goes, giving it a radiant and toxin-free appearance.



# Ultrasound

This piece is used to perform sonophoresis, dilating the pores through a cellular massage and then applying active products such as anti-acne peeling or moisturizers.



After the entire treatment, you will experience a feeling of relaxation. Ultrasound can release pain and unclog pores to treat fat.

# Advanced Cleaning

The newest way to hydrate your skin

Hydra-sh technology offers a unique and innovative skin rejuvenation treatment:

1. Skin cleansing, in which it is prepared for treatment.
2. Peeling, in which the active ingredient is placed on the skin and the pores are opened.
3. Suction pump or hydrodermabrasion, in which the superficial layer of the skin is removed (dead cells), impurities and blackheads are eliminated.
4. Sonophoresis, using ultrasound at 1.2 MHz, the anti-acne, anti-grease or anti-aging peel is introduced as required.

Hydra-sh accelerates the cell renewal process, resulting in smoother skin.

**Hydra-SH**



# Clinical Cases

Learn about some of the results obtained with Hydra-sh:



## **BLACKHEADS**

Courtesy: Tentrek Academy



## **COMEDONES**

Courtesy: Dr. Julieth Martinez



**BLACKHEADS**

Courtesy: Dr. Johvanna Tamayo



**ADVANCED CLEANING**

Courtesy: Dr. Maria Paula Gaitan



**BLACKHEADS**

Courtesy: Dr. Freddy Hernandez



**Hydra-SH**

# — Technical data sheet

<b>Ultrasound piece</b>	1.2 MHz
<b>Hydrodermabrasion piece</b>	Face (regular), nose (narrow)
<b>Microdermabrasion piece</b>	#140 / #170
<b>Suction</b>	Until 600mmHg ( $\pm 10\%$ )
<b>Consumption</b>	80Watts
<b>Reservoir capacity</b>	1200cc
<b>Feeding</b>	110VAC / 50 / 60Hz
<b>Dimensions (al*an*prof)</b>	1000x350x350mm
<b>Weight</b>	30kg

# Hydra-SH



REVIVE  
POWERED BY  TENTREK

